

A-dec® 300 Pro Cart

Ideal for both general and specialized practices, the A-dec 300 Pro Cart is configurable for your unique needs. And with its smooth, dynamic height adjustment, you can easily move from sitting to standing for all-day comfort as you work.

A-dec 343pro Cart

Features & Benefits:

- Clinical product integration from electric motors to intraoral scanners
- Dynamic height adjustment changes to fit your working position
- Gateway compatibility enables the A-dec+ software platform
- Optional assistant's arm accommodates 2- and 4-handed work
- Optional tray holder provides an additional worksurface
- Optional air and water quick connects support tabletop accessories
- Optional electrical duplex outlet adds 120 VAC power for tabletop accessories

Seamlessly integrate a full range of clinical devices including Primescan Connect™ from Dentsply-Sirona. Get connected with the A-dec+ software platform. Add the optional assistant's arm to accommodate both 2- and 4-handed procedures. And support your essential tabletop accessories. The A-dec 300 Pro Cart makes it all easy—and all your choice.

To learn more about the A-dec 300 Pro Cart, **contact your local authorized dealer.**

