

A-dec® 500 Pro

Government Package

A-dec 500 Pro supports efficiency and workflow without compromising design, durability, or patient comfort. Intuitive features on the responsive Dynamic Screen 7 (DS7) touchscreen including a simplified interface displaying only the information needed for the procedure at hand, and visual icons reflecting chair preset positions, further ensure that transferring personnel can quickly get up to speed.

Features:

- The A-dec 500 dental chair supports a patient weight capacity of 500 pounds (227 kg).
- Centrally located ambidextrous delivery system handle and touchscreen allows switching between left- and right-handed operators.
- Intuitive, highly configurable touchscreen interface lets doctors effortlessly switch between general restorative and endodontic functions.
- Six movable positions on the Traditional or Continental® delivery system allow operator to tilt the handpieces downward to easily grip when in use and return upright for safety.
- Customizable control head placement integrates most clinical devices.
- Drip-control vaporizer stops the last handpiece water droplet from falling onto the patient.
- Balanced flexarm with air brake and positive positioning moves smoothly and stays put.
- Signature A-dec ultra-thin back allows dentists to work in a healthy ergonomic position.
- Four adjustable points of rotation on the headrest support patient while providing an optimal view of the oral cavity from any angle.
- Sturdy multi-position armrests solidly support patients during entry and exit.
- When the chair reclines, the toeboard imperceptibly raises a full 9" (229 mm) for a gentle cradling effect.
- Water coolant indicator lights correspond with the handpiece in use for intuitive adjustments.
- Scientific mapping reduces pressure points and supports your patient's entire body.
- Precision hydraulic motion allows quiet, gentle stops and starts and eliminates jarring movements.
- Large 2-liter water bottle requires fewer refills.
- Supports 20+ user profiles for customized settings.
- Proudly manufactured in the USA and backed with a 10-year warranty.
- The A-dec+ connectivity option (Wi-Fi enabled) has been removed in adherence to federal government compliance requirements.



The A-dec 500 package includes:

- 511 dental chair with auto chair-stop
- 300-watt power supply
- 532/533 Pro delivery system with A-dec 500 DS7 touchscreen
- Programmable footswitch
- Removable tray holder
- User training and technical repair training/certification

Customize your A-dec 500 Pro package



Upholstery Sleek vacuum-formed upholstery eliminates seams for easy cleaning and maintenance.




Touchscreen Shown on the A-dec 532 Pro, the A-dec 500 DS7 touchscreen offers consistent and predictable intelligent control with a familiar touchscreen feel.




Delivery System A-dec 533 Pro offers precise positioning, corresponding LED coolant light indicators, and a configurable interface to switch between general restorative and endo functions.



Monitor Mounts Tilting mounts are VESA compliant. Available in Radius (shown), light post, wall, cabinet, ceiling or track mounted.



Stools A-dec 500 doctor's and assistant's stools offer a dynamic seating system that conforms and moves with the body. A-dec 400 line is also available.



Assistant's Holders A-dec 551sta provides built-in comfort and efficiency for assistants.



Assistant's Delivery A-dec 545sta assistant's instrumentation includes a self-contained 2-liter water bottle, 22.5" worksurface and multi-positionable instrumentation.



Dental Light 500 LED offers exceptional illumination and high color rendering index (CRI). Available in light post, wall, cabinet, ceiling or track mount.



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132
United States

Tel: 1.800.547.1883 Within USA/CAN
Tel: 1.503.538.7478 Outside USA/CAN
Fax: 1.503.538.0276
a-dec.com



The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2025 A-dec Inc. All rights reserved. 85.0836.00/4-25/Rev G