

# Ready for a next level delivery system? Save 35% when you upgrade your touchscreen!



Older A-dec 300 control pads

New A-dec Dynamic Screen 5 (DS5) touchscreen



Upgraded DS5 touchscreens show: chair position, RPM and torque, and apex locator (with Midwest)

[a-dec.com/promotions](https://a-dec.com/promotions)

Now through October 29, 2026, save 35% on a DS5 touchscreen with upgrade kit and unlock advanced functionality for your A-dec 300. Enjoy intuitive chair position icons, on screen messaging, multi-user settings, and command confirmation that helps prevent unintended actions—all designed to streamline workflow and elevate your clinical experience.

Talk to your A-dec representative about bringing next level control and confidence to your A-dec 300.

Part Number	Description	Promo Code
90.2204.00	Side Mount - with WiFi <b>Models:</b> 332, 333, 334, 335, 336, 342, 332pro, 333pro, 334pro, 335pro, 336pro, and 342pro	PR0239
90.2205.00	Post Mount - with Wi-Fi <b>Model:</b> 341pro	PR0239
90.2218.00	Post Mount - with Wi-Fi <b>Model:</b> 343pro	PR0239
90.2212.00	Wi-Fi Simulator 3 POSN <b>Models:</b> 41pro, 42pro	PR0242
90.2213.00	Wi-Fi Simulator 4 POSN <b>Models:</b> 41pro, 42pro	PR0242



## Compatibility Matrix

The chart below is meant to help you understand what additional upgrades are required for your delivery system to receive the full benefits of the DS5 and A-dec+

System Checks		Yes	No
<b>Delivery System</b>	<b>Is the delivery system a Continental® 300 non-300 Pro delivery system?</b>	Update the delivery to the pro-style whips.	No additional delivery conversions are required.
<b>Dental Chairs</b>	<b>Does the dental chair have a universal power supply and chair board? NOTE:</b> 311B/411B/511B dental chairs came standard with Universal Chair Board and Universal Power Supply after November 2023 (invoice should say “Pro compatible”). All chairs manufactured before this date will require additional upgrade kits for full DS5 functionality.	No additional chair conversions are required for full functionality.	In order to have utilization of all chair-related software features, the chair must be updated to a Universal Chair Board (UCB) and Universal Power Supply (UPS). Users who do not have a UCB may experience slower updates and the chair preset functions will not follow a multi-user account from room to room. Additional enhancements that involve chair functionality such as the occupancy sensors may not be available. <b>NOTE:</b> The 311A dental chair is not capable of upgrading to a UCB and UPS.

Your A-dec Customer Service team is available to answer your questions at 1.800.547.1883.

## What to Order

If additional upgrades are required, please find the appropriate part numbers and suggested retail price for the kits below.

Note: Additional upgrade kits are not covered under this promotional offer.

Updates	Order the Required Part Numbers:	Qty	Suggested Retail
300 Pro whip	90.1842.00 Whip upgrade	1	\$2,100
	Syringe tubing - <b>Choose one:</b>		
	23.1404.00 - Standard	1	\$57
	23.1301.01 - Warm water	1	\$765
	Handpiece tubing - <b>Select tubings:</b>		
	98.1162.00 - 4-hole	As needed	\$170
	98.1164.00 - 6-pin	As needed	\$245
	98.1163.00 - 5-hole fiber-optic with bulb	As needed	\$400
UCB and UPS update (Not available for 311A)	300w Universal Power Supply - <b>Choose one:</b>		
	90.1942.00 - 110-120V	1	\$990
	90.1943.00 - 220-240v	1	\$990
	Universal chair board		
	90.1963.00	1	\$1,070
	Digital inclinometers		
	90.1961.00	2	\$180 each
	511A chairs with pull-style potentiometer ONLY:		
	90.1639.00 - Inclinometer retrofit kit. NOTE: The installer will need to attach the new digital inclinometers to the bracket provided in 90.1639.00	1	\$440
	Wiring harness and hardware - <b>Choose one:</b>		
90.1962.00 - 311B/411	1	\$175	
90.1964.00 - 511A/511B	1	\$175	

Promotion applies only to the DS5 touchscreen upgrade kit. Depending on configuration, additional parts may be needed. Available in the U.S. and Canada.

\*Terms and Conditions. Offer is valid April 30–October 29, 2026, and must ship by December 31. All prices shown are MSRP and in U.S. dollars. Promo not available with new equipment purchase. Returns: Credit will be issued on unopened packages, minus 20% restocking fee and/or amount of qualified rebate. Promotional items are subject to product availability. A-dec reserves the right to change or cancel this offer at any time. This offer is void where prohibited, restricted, or taxed.