

A-dec® Choices
Which solution is right for you?

Choices



a dec®

There's an A-dec® for everyone.

A-dec offers reliable, creative solutions for practices of every size and specialty. Including yours. No matter which dental chair and delivery system you choose, you'll have the peace of mind of knowing they're engineered with legendary quality and reliability, and backed by a trusted A-dec warranty. So you can practise with confidence.



See a-dec.com/legal/warranty for details and exclusions.



A-dec 500

Premium A-dec performance, connectivity, and integration. Expand your practice possibilities while your patients enjoy superior comfort.

A-dec 500 Chair

- Exquisitely smooth and quiet, cradles patients with whole body support
- Impressive range of motion, including our lowest base-down position
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for optimal access and healthy, natural posture

A-dec 500 Pro Delivery

- Six adjustable control head positions for customisable placement
- DS7 touchscreen offers a simple, dynamic interface for intuitive, connected control
- Multi-user settings for personalised presets that follow you from operator to operator
- Supports implantology and advanced endodontics with integrated reciprocation, apex location, and file libraries*

* Availability varies by market



A-dec 400

A sleek design that brings exceptional patient comfort and dental team access together with consistent A-dec performance.

A-dec 400 Chair

- Distinctive contemporary style with contoured patient comfort
- Smooth, efficient positioning, plus 30° rotation on each side of centre
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for excellent access and healthy, natural posture

Delivery Options

- Pairs with A-dec 300 Pro delivery
- Choose from Continental® or Traditional styles
- Radius® design rotates around the chair—with no tools required—to support right- and left-handed clinicians



A-dec 300

Reliable, accessible and flexible, with A-dec+ connectivity to bring you full-size performance with a smaller footprint.

A-dec 300 Chair

- Compact, space-saving design fits a variety of operatory spaces
- Contoured, pressure-mapped surfaces support patients in comfort
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for excellent access and healthy, natural posture

A-dec 300 Pro Delivery

- Open platform design is easy to configure, customise, and upgrade in the future
- Pivoting control head reduces hand strain and fatigue
- Intuitive, responsive user interface shows just what you need to know in the moment
- Supports advanced endodontics with integrated reciprocation, apex location, and file libraries*



A-dec Performer eVo™

A compact, evolved solution with essential features for today's practices—all at a great value.

A-dec Performer eVo Chair

- Hydraulic lift ensures smooth, quiet, quick positioning
- Contoured surfaces and dual-articulating headrest for patient comfort
- Modern baseplate for easy installation and improved patient entry/exit
- Ergonomically designed for ample legroom and excellent clinician access

A-dec 300 Delivery

- Four handpiece holder positions allow for customisable configurations
- Balanced flex arm with air brake for easy positioning and stability
- Choice of Traditional or Continental® style
- Radius® design rotates to right or left side of chair, with no tools required

Chair**	A-dec 500	A-dec 400	A-dec 300	A-dec Performer eVo
Ultra-soft start/stop hydraulic drive	•	—	—	—
Hybrid drive (hydraulic base/electric tilt)	—	•	•	•
Swivel (30° each side of center)	•	•	—	•
Virtual pivot for added patient comfort	•	•	•	•
Thin backrest 1.25" (32mm)	•	•	•	•
Dual-articulating gliding headrest	Lever release	Locking knob or lever release	Locking knob	Locking knob
Multi-position armrest settings	2	2	2	2
Chair at lowest point	343 mm (13.5")	349 mm (13.75")	349 mm (13.75")	398 mm (15.7")
Chair at highest point	800 mm (31.5")*	800 mm (31.5")	800 mm (31.5")	835 mm (32.9")
Integrated 300 watt power supply	•	•	•	•
Load rating	500 lb. (227 kg)	400 lb. (181 kg)	400 lb. (181 kg)	400 lb. (181 kg)
Choice of sewn or formed upholstery	•	•	•	Formed (available in select colours)
Dental light	A-dec 500 LED	A-dec 500 LED, A-dec 300 LED	A-dec 500 LED, A-dec 300 LED	Performer eVo LED

Delivery System**	A-dec 500 Pro	A-dec 300 Pro	A-dec 300
Chairs	A-dec 500	A-dec 500, A-dec 400, A-dec 300	A-dec 500, A-dec 400, A-dec 300, Performer eVo
Mounting options	Radius	Radius, side, and support centre	Radius, side, and support centre
Control pad	DS7	CP5, CP5i, DS5	Deluxe, Standard, or no touchpad
A-dec+ connectivity	•	Optional (requires A-dec+ Gateway or DS5)	—
Four-position control block	•	Optional	Optional
Three-position control block	—	•	•
Holder positions	6	5	5
Clinical device integration	5	2	2
USB pass-through port	•	(on deliveries with CP5i)	—
Support centre	A-dec 400	A-dec 400, A-dec 300	A-dec 400, A-dec 300
Self-contained water system with quick-disconnect 2-litre bottle	•	•	•
Accommodates amalgam separator	•	•	•
Cuspidor swivel	•	A-dec 400 only	A-dec 400 only
Programmable timing for cup and bowl	•	•	•

*Using maximum upper limit. The factory set upper limit is 737 mm (29").

** Availability varies by market; check with your A-dec Territory Manager or authorised A-dec dealer.

Reliable Creative Solutions

FOR DENTAL PROFESSIONALS WORLDWIDE




MORE THAN A WARRANTY. IT'S A COMMITMENT.

At A-dec, our focus is on making dentistry better—which includes superior customer service along with superior dental products. Your A-dec operatory equipment is backed by a ten-year warranty, and a lifetime of customer support. We'll be here today, tomorrow, and well into the future. That's our commitment to you.

For complete details, including exclusions, visit: a-dec.com/legal/warranty

Any colours and products shown may vary and are subject to change without notice.
Contact your authorised A-dec dealer for the most current product information.
All trademarks in this document are subject to the trademark rights of A-dec, its affiliates, and their collective suppliers in the United States and foreign countries.

©2026 A-dec Inc. All rights reserved. 85.1093.25/BB/1M/3-26/Rev B



A-dec Headquarters

2601 Crestview Drive

Newberg, Oregon 97132 USA

P: 1.800.547.1883 within USA/CAN

P: +1.503.538.7478 outside USA/CAN

a-dec.com

A-dec Australia

P: 1.800.225.010 within AUS

P: +61.(0)2.8332.4000 outside AUS

A-dec United Kingdom

P: 0800.233285 within UK

P: +44.(0).24.7635.0901 outside UK

