

# Performer eVo™

Legendary A-dec® reliability and performance  
come together in a smart, evolved solution.



a dec®

# A-DEC® QUALITY AT AN IMPRESSIVE VALUE.

Performer eVo™ brings legendary A-dec reliability and performance together in a smart solution made for the way you practise. Get the essential features you want, the comfort your patients deserve, and the dependable quality A-dec is known for.

Designed for minimal maintenance, and backed by a trusted A-dec warranty, Performer eVo delivers lasting value day after day, patient after patient. It's the best of A-dec, thoughtfully evolved.

## Optional upgrades:

- Deluxe touchpad
- Monitor mount
- Doctor's and Assistant's stools
- Cuspidor and support centre
- Assistant's instrumentation
- Sewn upholstery
- Clinical devices



See [a-dec.com/legal/warranty](https://www.a-dec.com/legal/warranty) for details and exclusions.

Any colours and products shown may vary and are subject to change without notice. Contact your authorised A-dec dealer for the most current product information. All trademarks in this document are subject to the trademark rights of A-dec, its affiliates, and their collective suppliers in the United States and foreign countries.

Not available in all countries. Check for local availability.  
©2026 A-dec Inc. All rights reserved. 85.1090.25/BB/500/2-26/Rev A

## Performer eVo chair features:

- Compact design—fits a wide variety of operatory spaces
- 60° chair swivel with brake—for easy rotation and patient positioning
- Dual-articulating headrest—comfortably supports patient's head and neck
- Hydraulic lift—for smooth, quiet, quick positioning
- Modern baseplate—for easy installation and improved patient entry/exit.
- Patient weight limit 181kg /400 lbs

## A-dec 300 delivery system features:

- Choice of Traditional or Continental® style
- Radius® mounted—easily rotates around chair without tools, for right- and left-handed compatibility
- Balanced flex arm with air brake—for easy positioning and stability
- Four handpiece holder positions—customisable configurations
- 2-litre, quick disconnect, self-contained water bottle—to easily control treatment water quality
- Adjustable air coolant—prevents handpiece damage from frictional heat
- Individually adjustable water coolant—accommodates coolant needs for each handpiece

## Performer eVo LED light features:

- Low heat/no noise—for patient and doctor's comfort
- Three light intensities—to support a range of procedures
- LED technology—for accurate diagnosis, reduced eyestrain and unparalleled ergonomics
- Touchless on/off sensor—easier to clean with fewer barriers needed
- Cure-safe mode—prevents premature curing of dental composites
- Three-axis head—for optimal positioning

## The Performer eVo package:



Performer eVo chair



A-dec 300 delivery system



Performer eVo LED light



Standard touchpad




Performer eVo doctor's stool

## Choose from three formed upholstery colours:



Sandstone



Crow



Olympic