



Porta-Chair

A transportable, quick-setup dental solution.

Even temporary setups require full functionality. Porta-Chair delivers.


A-dec® has been designing portable dental chairs to bring oral care to underserved communities and remote locations since 1974. With its unique flip-down toeboard design, adjustable seat height, and chair back angle, the updated Porta-Chair provides optimal oral cavity access—no matter where you are.

A small chair with big benefits.

A built-in toeboard with knee-break functionality, faster set up time, and 350-pound patient weight capacity allow you to serve patients professionally and comfortably. When you've completed your clinic, the entire chair neatly folds into its hard-shell carrying case for easy transport.

Chair features:

- **48 POUNDS** — Lightweight enough for one person to transport.
- **350 POUND PATIENT WEIGHT CAPACITY** — Accommodates patients from pediatric to adult.
- **ADJUSTABLE SEAT** — Adjustable from 21" to 26" for patient and caregiver comfort.
- **ADJUSTABLE CHAIR BACK** — Angles from 0° to 80° adapts to different shapes and sizes.
- **IMPROVED CARRYING CASE** — Durable protective case with wheels for easy transportation.
- **FULLY INTEGRATED TOE BOARD** — Easily flips down and folds without additional components.
- **KNEE-BREAK FEATURE** — Toeboard can be dropped into knee-break position to aid chair entry and exit.
- **DURABLE UPHOLSTERY** — Soft, formed upholstery is durable and easy to clean.



A-dec PAC 1 Portable Dental System

Complete your portable treatment capability with the field-proven unit that offers manual control for two handpieces, clip on water and vacuum systems, and a self-contained water system.



Contact customer.service@a-dec.com or 1.800.547.1883 for more information.