

Cavitron.

For dental professionals searching for the  
streamlined look of integration,

**Cavitron® Integrated Systems  
gives you the freedom of  
more open space.**

### Cavitron® Integrated Systems

- Fewer cords and less visible equipment
- Open office environment preserves counter space
- Compatible with Steri-Mate® 360° Rotating Handpiece



**There's only one Cavitron®.**

**Cavitron®**  
**integrated**  
**systems**

## Cavitron® Integrated Systems offer an open environment with the same benefits of traditional Cavitron® units.

- Compatible with **Steri-Mate® 360°**, a rotating handpiece that allows the clinician to fully rotate an ultrasonic insert without stopping for adjustments
- Detachable and autoclavable handpiece reduces cross contamination for enhanced infection control
- **Patented Sustained Performance System**, SPS™, **Technology** automatically adjusts power to maintain scaling efficiency, even at low power
- **BlueZone™** extended low power range designed for improved patient comfort



## Cavitron Integrated Systems integrate easily with A-dec® 500 and A-dec® 300 delivery systems.

For dental professionals searching for the streamlined look of integration, you are now able to combine Cavitron Systems with the convenience of integration into the A-dec 500 or A-dec 300 dental units.

| Item                                          | A-dec Integrated Part # | Description                       | Compatible With:                             |
|-----------------------------------------------|-------------------------|-----------------------------------|----------------------------------------------|
| Cavitron A-dec Built-In Ultrasonic Scaler Kit | 54.0599.00              | Cavitron Steri-Mate 360° Kit, 61" | A-dec 532, 533, 542, 332, 333, 334, 335, 336 |
|                                               | 54.0597.00              | Cavitron Steri-Mate Kit, 61"      |                                              |
|                                               | 54.0600.00              | Cavitron Steri-Mate 360° Kit, 82" | A-dec 541                                    |
|                                               | 54.0598.00              | Cavitron Steri-Mate Kit, 82"      |                                              |

### For an integrated Cavitron kit contact:

1.800.547.1883 Within USA  
1.503.538.7478 Outside USA/Canada  
www.a-dec.com



### For a built-in Cavitron kit for field installation contact:

Dentsply Sirona  
1.800.989.8826  
dentsplysirona.com

